

**Pronto®
Configurations Quick
Reference Guide**

Designed to help you select the right Pronto model for your client. All Pronto wheelchairs feature SureStep® Technology.

Pronto® Series



Group 2

Group 3



M41™

The Pronto M41 features MK₅™ electronics and high efficiency motors for longer battery range and a faster top speed. The M41 is available with either a lightweight fold-down seat or a semi-recline Captain's seat with width and height adjustable arms. The M41 also features a 300 lb. weight capacity, 5 mph top speed, 10.5" x 3.5" drive wheels and U1 batteries.



M51™/M51-R/M51CG

The Pronto M51 is designed to meet the needs of a wide range of consumers. The M51 model comes standard with a center-post Captain's seat and the M51-R comes standard with a four-post Rehab seat. The M51CG comes standard with a center-of-gravity power tilt system. All models have a standard 300 lb. weight capacity, 4.25 mph top speed, 10.5" X 3.5" drive wheels, and U1 batteries.



M61™/M61-R

The Pronto M61 comes standard with a Power Elevating Seat that offers up to 5" of lift with the flick of a joystick. The M61 model comes standard with a Captain's seat and the M61-R comes standard with a Rehab seat. The chair also features a unique building block system that allows for different starting seat-to-floor heights. Both models have a standard 300 lb. weight capacity, 4.25 mph top speed, 10.5" X 3.5" drive wheels and U1 batteries.



M91™/M91-R/M91-TS

The M91-TS for the Group 3 single power option code features a Power Tilt System (PTO Plus) and it offers upgradeable electronics, alternative driver controls, and switch options. It is available in both a 300 lb. weight capacity (top speed of 6 mph) and a 400 lb. weight capacity with the Heavy Duty package (top speed 5 mph). The M91 model comes standard with a Captain's seat and the M91-R comes standard with an Adjustable Rehab seat. All models come standard with Traction Control Design (TCD), Stability Lock, 14" X 3" drive wheels and 22NF batteries.



M94™

The Pronto M94 has a weight capacity of 500 lb. offering a complete Bariatric Solution while maintaining a small wheelbase and a 19.5" turning radius. Seat sizes range from 18" to 24" wide on both the Captain's seat and Rehab seat to accommodate a wide range of consumers. The M94 comes standard with TCD, 4.5 mph top speed, Stability Lock, 14" X 3" drive wheels, electronics and 22NF batteries.

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E., Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2008 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Form No. 05-134 Rev. 05/08

	M41™	M51™/M51-R /M51CG	M61™/M61-R	M91™/M91-R/M91-TS	M94™
Weight Capacity	300 lb.	300 lb.	300 lb.	M91: 300 lb.; M91-R: 400 lb. M91-TS:275 lb.; M91-TS: 375 lb. w/powerd seating	500 lb.
Speed	5 mph	4.25 mph	4.25 mph	6 mph;5 mph (M91-R, M91-TS)	5 mph
Motor Package	2-pole in-line	2-pole	2-pole	4-pole/4-pole Heavy Duty	4-pole Heavy Duty
Seat Width	16", 18", 20"	16", 18", 20"	16", 18", 20"	16", 18", 20", 22"	18", 20", 22", 24"
Seat Depth	16" - 18"	16" - 20"	16" - 20"	16" - 22"	16" - 22"
Electronics	MK ₅ ™ NX™ with SPJ™+ Joystick	MK ₅ NX with SPJ+ Joystick	M61: MK ₅ NX w/acc M61-R: MK ₅ SPJ + w/acc MK660-TS: w/acc.SPJ + w/acc	M91/M91 HD: MK ₅ NX-75 w/SPJ+ M91R/M91-R HD: MK690 w/SPJ+ M91-TS/M91-TS HD: MK690 ACC w/SPJ+ w/PSS	MK690 / SPJ+
Tire/Wheel Size	10.5" x 3.5"	10.5" x 3.5"	10.5" x 3.5"	14" x 3"	14" x 3"
Batteries	U1	U1	U1	22NF	22NF
Turning Radius	17.75" Front - 19.5" Rear	19.5"	19.5"	21.5"	21.5"
Base Width	24"	24"	24"	25.9"	25.9"
Base Length	34"	34"	34"	39"	39"
STF Height*					
Low	18"	17"	17"	18.5"	18.5"
Medium	20"	19"	19"	21"	21"
Tall	22"	21"	21"	23"	23"
Core Technologies					
Center-Wheel Drive	Standard	Standard	Standard	Standard	Standard
SureStep®	Standard	Standard	Standard	Standard	Standard
Stability Lock	NA	NA	NA	Standard	Standard
MK₅™ Electronics	Standard	Standard	Standard	Standard	Standard
Traction Control Design (TCD)	NA	NA	NA	Standard	Standard
Transportability**	NA	Transport Ready Option	NA	Transport Ready Option	NA
Powered Seating Options	NA	CG Tilt	Standard Power Elevating Seat	Formula PTO Plus	NA
SADMERC Group Level	2	2	2 Seat Elevator	3 / 3HD / 3 Single Power / 3HD Single Power	3 VHD
SADMERC Placement Code by Model	K0823	K0823 - M51 K0822 - M51-R K0835 - M51-CG	K0831 - M61 K0830 - M61-R	K0849 - M91 K0848 - M91- R K0851 - M91- HD K0856 - M91- TS K0850 - M91- R HD K0858 - M91- HD K0858 - M91- TS HD	K0853 - M94 K0852 - M94-R

* Seat-to-floor heights are approximate and vary dependent upon seat depth.

** Invacare recommends that wheelchair users not be transported in vehicles of any kind while in wheelchairs. As of this date, the U.S. Department of Transportation has not approved any tie-down system for transportation of a user while in a wheelchair in a moving vehicle of any type.

NOTICE: This information is not intended to be, nor should be, considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor when needed to discuss specific situations in further detail.

Information current as of April 2008